

**Safety Data Sheet**  
acc. to OSHA HCS

Printing date 12/15/2022

Reviewed on 12/15/2022

**1 Identification**

- **Product identifier**
- **Trade name:** Opalescence Go™ 10% (Mint and Melon)
- **Article number:** SDS 481-001.01, 1005956, 1005957
- **Application of the substance / the mixture** Professional Dental Teeth Bleaching Gel
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
Ultradent Products Inc.  
505 W. Ultradent Drive (10200 S)  
South Jordan, UT 84095-3942  
USA  
onlineordersupport@ultradent.com
- **Information department:** Customer Service
- **Emergency telephone number:**  
CHEMTREC (NORTH AMERICA) : (800) 424-9300  
(INTERNATIONAL) : +(703) 527-3887

**2 Hazard(s) identification**

- **Classification of the substance or mixture**



Sensitization - Skin 1 H317 May cause an allergic skin reaction.

- **Label elements**
- **GHS label elements** The product is classified according to the Globally Harmonized System (GHS)
- **Hazard pictograms** GHS07
- **Signal word** Warning
- **Health Hazard-determining components of labeling:**  
Artificial Watermelon  
Oils, Peppermint
- **Hazard statements**  
H317 May cause an allergic skin reaction.
- **Precautionary statements**  
P261 Avoid breathing dust/fume/gas/mist/vapors/spray  
P272 Contaminated work clothing must not be allowed out of the workplace.  
P280 Wear protective gloves.  
P302+P352 If on skin: Wash with plenty of water.  
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.  
P321 Specific treatment (see on this label).  
P363 Wash contaminated clothing before reuse.  
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

(Contd. on page 2)

# Safety Data Sheet

acc. to OSHA HCS

Printing date 12/15/2022

Reviewed on 12/15/2022

**Trade name: Opalescence Go™ 10% (Mint and Melon)**

(Contd. of page 1)

· **Classification system:**

· **NFPA ratings (scale 0 - 4)**



· **HMIS-ratings (scale 0 - 4)**

HEALTH	1	Health = 1
FIRE	0	Fire = 0
REACTIVITY	0	Reactivity = 0

### 3 Composition/information on ingredients

· **Chemical characterization: Mixtures**

· **Description:** Mixture of the substances listed below with nonhazardous additions.

· **Dangerous components:**

56-81-5	Glycerin	>30-<50%
7722-84-1	Hydrogen Peroxide	10%
	Potassium Hydroxide	>1-<5%
	Artificial Watermelon	>1-<5%
	Silicon Dioxide Chemically Prepared	>1-<5%
	Dipotassium Phosphate	>1-<5%
8006-90-4	Oils, Peppermint	<1%
7681-49-4	Sodium Fluoride	<1%

· **Additional information:**

The specific chemical identity of composition is being withheld as a trade secret. The specific chemical identity is made available to health professionals, employees, and designated representatives in accordance with the applicable provisions of paragraph §1910.1200.

### 4 First-aid measures

· **Description of first aid measures**

· **After inhalation:**

Supply fresh air and to be sure call for a doctor.

In case of unconsciousness place patient stably in side position for transportation.

· **After skin contact:** Immediately wash with water and soap and rinse thoroughly.

· **After eye contact:** Rinse opened eye for several minutes under running water.

· **After swallowing:** If symptoms persist consult doctor.

· **Information for doctor:**

· **Most important symptoms and effects, both acute and delayed** No further relevant information available.

· **Indication of any immediate medical attention and special treatment needed**

No further relevant information available.

### 5 Fire-fighting measures

· **Extinguishing media**

· **Suitable extinguishing agents:** Use fire fighting measures that suit the environment.

(Contd. on page 3)

# Safety Data Sheet

acc. to OSHA HCS

Printing date 12/15/2022

Reviewed on 12/15/2022

**Trade name: Opalescence Go™ 10% (Mint and Melon)**

(Contd. of page 2)

- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** No special measures required.

## 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:**
  - Dilute with plenty of water.
  - Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**
  - Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
  - Dispose contaminated material as waste according to item 13.
  - Ensure adequate ventilation.
- **Reference to other sections**
  - See Section 7 for information on safe handling.
  - See Section 8 for information on personal protection equipment.
  - See Section 13 for disposal information.

## 7 Handling and storage

- **Handling:**
- **Precautions for safe handling**
  - See product labeling.
  - Ensure good ventilation/exhaustion at the workplace.
  - Prevent formation of aerosols.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
  - Requirements to be met by storerooms and receptacles:** No special requirements.
  - Information about storage in one common storage facility:** Not required.
  - Further information about storage conditions:** None.
  - Specific end use(s)** No further relevant information available.

## 8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**
  - The following constituents are the only constituents of the product which have a PEL, TLV or other recommended exposure limit.
  - At this time, the other constituents have no known exposure limits.

### 56-81-5 Glycerin

PEL	Long-term value: 15* 5** mg/m <sup>3</sup> mist; *total dust **respirable fraction
TLV	TLV withdrawn-insufficient data human occup. exp.

### 7722-84-1 Hydrogen Peroxide

PEL	Long-term value: 1.4 mg/m <sup>3</sup> , 1 ppm
-----	--

(Contd. on page 4)

# Safety Data Sheet

acc. to OSHA HCS

Printing date 12/15/2022

Reviewed on 12/15/2022

**Trade name: Opalescence Go™ 10% (Mint and Melon)**

(Contd. of page 3)

REL	Long-term value: 1.4 mg/m <sup>3</sup> , 1 ppm
TLV	Long-term value: 1 ppm A3
<b>Potassium Hydroxide</b>	
REL	Ceiling limit value: 2 mg/m <sup>3</sup>
TLV	Ceiling limit value: 2 mg/m <sup>3</sup>
<b>Silicon Dioxide Chemically Prepared</b>	
TWA	Short-term value: 0.8 mg/m <sup>3</sup>

· **Additional information:** The lists that were valid during the creation were used as basis.

· **Exposure controls**

· **Personal protective equipment:**

· **General protective and hygienic measures:**

Immediately remove all soiled and contaminated clothing.

Wash hands before breaks and at the end of work.

· **Breathing equipment:**

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

· **Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material is based on consideration of the penetration times, rates of diffusion and the degradation

· **Material of gloves**

The selection of suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye protection:** Goggles recommended during refilling.

· **Body protection:** Protective work clothing

## 9 Physical and chemical properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

Form: Gel

Color: White Opaque

· **Odor:** Mint

Melon

· **Odor threshold:** Not determined.

· **pH-value at 20 °C:** 5-7

(Contd. on page 5)

# Safety Data Sheet

acc. to OSHA HCS

Printing date 12/15/2022

Reviewed on 12/15/2022

Trade name: Opalescence Go™ 10% (Mint and Melon)

(Contd. of page 4)

· <b>Change in condition</b>	
Melting point/Melting range:	Undetermined.
Boiling point/Boiling range:	Undetermined
· <b>Flash point:</b>	Not applicable
· <b>Flammability (solid, gaseous):</b>	Not applicable.
· <b>Decomposition temperature:</b>	Not determined.
· <b>Auto igniting:</b>	Product is not selfigniting.
· <b>Danger of explosion:</b>	Product does not present an explosion hazard.
· <b>Explosion limits:</b>	
Lower:	Not determined.
Upper:	Not determined.
· <b>Vapor pressure:</b>	Not determined.
· <b>Density at 20 °C:</b>	1.25 g/cm <sup>3</sup>
· <b>Relative density</b>	Not determined
· <b>Vapor density</b>	Not determined.
· <b>Evaporation rate</b>	Not determined.
· <b>Solubility in / Miscibility with Water:</b>	Partly miscible.
· <b>Partition coefficient (n-octanol/water):</b>	Not determined.
· <b>Viscosity:</b>	
Dynamic:	Not determined.
Kinematic:	Not determined
· <b>Other information</b>	No further relevant information available.

## 10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

## 11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**

· **LD/LC50 values that are relevant for classification:**

**ATE (Acute Toxicity Estimate)**

Oral	LD50	2,684 mg/kg
Dermal	LD50	70,000 mg/kg (rat)
Inhalative	LC50/4 h	>39.7 mg/l

(Contd. on page 6)

US

# Safety Data Sheet

acc. to OSHA HCS

Printing date 12/15/2022

Reviewed on 12/15/2022

**Trade name: Opalescence Go™ 10% (Mint and Melon)**

(Contd. of page 5)

<b>56-81-5 Glycerin</b>		
Oral	LD50	7,750 mg/kg (Guinea pig) 4,100 mg/kg (mouse) 5,570 mg/kg (rat) 27,000 mg/kg (rabbit)
Dermal	LC50 Fish	>5,000 mg/l (Fish)
	LD50	>21,900 mg/kg (rat) 10,000 mg/kg (rabbit)
<b>7722-84-1 Hydrogen Peroxide</b>		
Oral	LC50 Fish	16.4 mg/l (Fish)
<b>Potassium Hydroxide</b>		
Oral	LD50	214 mg/kg (rat)
	LC50 Fish	80 mg/l (Fish)
<b>Silicon Dioxide Chemically Prepared</b>		
Oral	LD50	>15,000 mg/kg (mouse) >3,300 mg/kg (rat)
	LC50 Fish	>10,000 mg/l (Fish) (Toxicity to fish)
Dermal	LD50	>5,000 mg/kg (rabbit)
Inhalative	LC50/4 h	0.139 mg/l (rat)
<b>Dipotassium Phosphate</b>		
Oral	LD50	4,260-5,700 mg/kg (rat)
Dermal	LD50	>5,000 mg/kg (rabbit)
Inhalative	LC50/4 h	>0.83 mg/l (rat)
<b>8006-90-4 Oils, Peppermint</b>		
Oral	LD50	2,490 mg/kg (mouse)
		2,426 mg/kg (rat)
<b>7681-49-4 Sodium Fluoride</b>		
Oral	LD50	52 mg/kg (mouse)
	LC50 Fish (static)	17 mg/l (Fish)
Dermal	LD50	175 mg/kg (rat)

· **Primary irritant effect:**

· **on the skin:** No irritant effect.

· **on the eye:** No irritating effect.

· **Sensitization:** Sensitization possible through skin contact.

· **Additional toxicological information:**

The product shows the following dangers according to internally approved calculation methods for preparations:  
Irritant

· **Carcinogenic categories**

· <b>IARC (International Agency for Research on Cancer)</b>		
7722-84-1	Hydrogen Peroxide	3
9003-39-8	Polyvinylpyrrolidone	3
7681-49-4	Sodium Fluoride	3
79-10-7	acrylic acid	3

(Contd. on page 7)

# Safety Data Sheet

acc. to OSHA HCS

Printing date 12/15/2022

Reviewed on 12/15/2022

**Trade name: Opalescence Go™ 10% (Mint and Melon)**

(Contd. of page 6)

· **NTP (National Toxicology Program)**

None of the ingredients is listed.

· **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.

## 12 Ecological information

· **Toxicity**

· **Aquatic toxicity:**

**56-81-5 Glycerin**

EC50	>10,000 mg/kg (Bacteria)
------	--------------------------

**7722-84-1 Hydrogen Peroxide**

EC50	1.38 mg/l (Algae) 2.4 mg/l (daphnia)
------	---

**Silicon Dioxide Chemically Prepared**

EC50	>1,000 mg/kg (daphnia)
------	------------------------

**7681-49-4 Sodium Fluoride**

EC50	272 mg/kg (Algae) 98 mg/kg (daphnia)
Algae Toxicity (static)	7 mg/l (Algae)

· **Persistence and degradability** No further relevant information available.

· **Behavior in environmental systems:**

· **Bioaccumulative potential** No further relevant information available.

· **Mobility in soil** No further relevant information available.

· **Additional ecological information:**

· **General notes:**

Water hazard class 1 (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

· **Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

· **Other adverse effects** No further relevant information available.

## 13 Disposal considerations

· **Waste treatment methods**

· **Recommendation:**

Dispose of contents/container in accordance with international, federal, state, and local regulations.

· **Uncleaned packagings:**

· **Recommendation:** Disposal must be made according to official regulations.

· **Recommended cleansing agent:** Water, if necessary with cleansing agents.

US

(Contd. on page 8)

# Safety Data Sheet

acc. to OSHA HCS

Printing date 12/15/2022

Reviewed on 12/15/2022

Trade name: Opalescence Go™ 10% (Mint and Melon)

(Contd. of page 7)

## 14 Transport information

<b>· UN-Number</b>	
<b>· DOT, IMDG, IATA</b>	Not Regulated
<b>· UN proper shipping name</b>	
<b>· DOT, IMDG, IATA</b>	Not Regulated
<b>· Transport hazard class(es)</b>	
<b>· DOT, ADN, IMDG, IATA</b>	
<b>· Class</b>	Not Regulated
<b>· Packing group</b>	
<b>· DOT, IMDG, IATA</b>	Not Regulated
<b>· Environmental hazards:</b>	Not Applicable.
<b>· Special precautions for user</b>	Not Applicable
<b>· Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code</b>	Not Applicable.
<b>· UN "Model Regulation":</b>	Not Regulated

## 15 Regulatory information

<b>· Safety, health and environmental regulations/legislation specific for the substance or mixture</b>		
No further relevant information available.		
<b>· Sara</b>		
<b>· Section 355 (extremely hazardous substances):</b>		
7722-84-1	Hydrogen Peroxide	
<b>· Section 313 (Specific toxic chemical listings):</b>		
7757-79-1	Potassium Nitrate	
<b>· TSCA (Toxic Substances Control Act):</b>		
56-81-5	Glycerin	ACTIVE
7722-84-1	Hydrogen Peroxide	ACTIVE
	Potassium Hydroxide	ACTIVE
	Dipotassium Phosphate	ACTIVE
8006-90-4	Oils, Peppermint	ACTIVE
7681-49-4	Sodium Fluoride	ACTIVE
<b>· Hazardous Air Pollutants</b>		
None of the ingredients is listed.		
<b>· Proposition 65</b>		
<b>· Chemicals known to cause cancer:</b>		
None of the ingredients is listed.		
<b>· Chemicals known to cause reproductive toxicity for females:</b>		
None of the ingredients is listed.		
<b>· Chemicals known to cause reproductive toxicity for males:</b>		
None of the ingredients is listed.		

(Contd. on page 9)



# Safety Data Sheet

acc. to OSHA HCS

Printing date 12/15/2022

Reviewed on 12/15/2022

**Trade name: Opalescence Go™ 10% (Mint and Melon)**

(Contd. of page 8)

· **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

None of the ingredients is listed.

· **ACGIH Carcinogenicity (American Conference of Governmental Industrial Hygienists)**

7722-84-1	Hydrogen Peroxide	A3
7681-49-4	Sodium Fluoride	A4

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **Chemical safety assessment:** A chemical safety assessment has not been carried out.

## 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:** Environmental, Health, and Safety

· **Contact:** Customer Service

· **Date of preparation / last revision** 12/15/2022 / -

· **Abbreviations and acronyms:**

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

Sensitization - Skin 1: Skin sensitisation – Category 1

· **\* Data compared to the previous version altered.**